



Major Repair Kit CAT. No. 1400K-04 for Model 1400

ID	Qty.	Description	ID	Qty.	Description
Major Repair Kit 1400K-04			12	1	Intake Screen
1	1	Mechanical Shaft Seal Kit, 5/8" CAT. No. 1401E	13	1	O-Ring
2	1	Gasket Oil Case	14	1	O-Ring
3	1	Lip Seal	15	2	Woodruff Key
4	1	Exhaust Valve	16	1	Retaining Ring
5	1	Stud for Valve Cover	17	1	Woodruff Key
6	1	Flat Washer Brass			
7	4	Vane Spring			
8	2	Spring Holder for Vane			
9	4	Vane			
10	1	O-Ring	Minor Repair Kit 1400K-03		
11	1	Gasket, Intake Chamber Cover	Contains same as Major Kit 1400K-04 except vane ID No. 9.		

DuoSeal 1400 Repair Kits