



DuoSeal 1400 Spare Parts

ID	Qty.	CAT. No.	Description	ID	Qty.	CAT. No.	Description
A	1	41-2191	Pump Pulley with Set Screw			★1417	Standard Exhaust Filter
B	1	61-8540	Oil Case Assembly			★1417L	Replacement Exhaust Filter Element (2 pk_
C	1	41-1045	Front Plate			★1407K-11	DuoSeal Vacuum Pump Oil -1 Qt.
D	1	41-1042	Inner Ring			★1407K-15	DuoSeal Vacuum Pump Oil -1 Ga.
E	1	41-1040	Rotor			★331040-5	5' x 7/16" I.D. Vacuum Hose
F	1	61-8702	Back End Plate	9	1	41-1059	Vane
G	1	1393J	Hose Adapter for 7/16" I.D. Hose		1	61-8901	Belt Guard
H	1	61-8675	Intake Ring		1	61-8916	Pump Base
I	1	41-1060	Shaft	18	1	41-1736	Gas Ballast Valve
					1	41-0551	Motor Pulley with Set Screw
1	1	★1401E	Mechanical Shaft Seal Kit, 5/8"	19	1	726130	Oil Drain Valve
	1	★1401E-01	Lip Shaft Seal Kit, 5/8"		1	61-8622	1/3 HP Motor, 115V, 60Hz, 1Ph
	1	★1400A	V-Belt				