



DuoSeal 1402B-01 Spare Parts

ID	Qty.	CAT. No.	Description
A	1	61-8101	Oil Case Assembly
B	1	1393K	Intake Adapter for 13/16" I.D. Hose
C	1	41-2045	Small End Plate
D	1	61-8108D	Exhaust Ring
E	1	41-1518	Exhaust Rotor
F	1	41-2074	Pump Pulley with Set Screw
G	1	41-2093	Shaft
H	1	61-8250C	Center Plate
I	1	41-2092	Intake Rotor
J	1	61-8113C	Intake Ring

ID	Qty.	CAT. No.	Description
1	1	★1401D	Mechanical Shaft Seal Kit, 3/4"
	1	★1405A	V-Belt
		★1417P-10	Standard Exhaust Filter
		★1417G	Replacement Exhaust Filter Element
		★1407K-11	DuoSeal Vacuum Pump Oil -1 Qt.
		★1407K-15	DuoSeal Vacuum Pump Oil -1 Ga.
		★331060-5	5' x 13/16" I.D. Vacuum Hose

ID	Qty.	CAT. No.	Description
18	1	41-1487	Vane, Small (Requires 2)
19	1	41-0959	Vane, Large (Requires 2)
	1	1371A-01	Belt Guard & Base Kit
	1	61-8803	Belt Guard/Cover
23	1	41-1736	Gas Ballast Valve
	1	41-0668	Motor Pulley with Set Screw
24	1	61-8912	Exhaust Plug/Cover
	1	1401D-01	Shaft Lip Seal Kit, 3/4"
	1	61-8698	Pump Base
	1	41-1907A	1/2 HP Motor, 115/230V, 60Hz, 1Ph
30	2	61-2229A	Oil Window Glass (Requires 2)
31	1	726130	Drain Valve
32	2	61-2370	Oil Window Cover (Requires 2)
33	4	61-1014	Screw (Requires 4)
34	2	66-0304	O-Ring (Requires 2)